

— ABLYBIO, Help Your Research



# Histone H3 (Mono Methyl Lys4) Rabbit mAb

货号: **AA Y1108**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA CHIP Cut&Tag
推荐浓度	<b>IHC:</b> 1:5000-1:20000 <b>WB:</b> 1:2000-1:10000 <b>IF:</b> 1:200-1:1000 <b>ELISA:</b> 1:5000-1:20000 <b>IP:</b> 1:50-1:200 <b>CHIP:</b> 1:50-1:100 <b>Cut&amp;Tag:</b> 1:50-1:100
理论分子量	15kD
实测分子量	15kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.110.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	
基因ID	
基因名	HIST1H3A/HIST1H3/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B/H3F3C;Histone H3.1/Histone H3.2/Histone H3.3/Histone H3.3 C;H3K27ME2;HIST1H3A;H3FA;HIST1H3B;H3FL;HIST1H3C;H3FC;HIST1H3D;H3FB;HIST1H3E;H3FD;HIST1H3 F;H3FI;HIST1H3G;H3FH;HIST1H3H;H3FK;HIST1H3I;H3FF;HIST1H3J;H3FJ;Histone H3.1;Histone H3/a;Histone H3/b;Histone H3/c;Histone H3/d;Histone H3/f;Histone H3/h;Histone H3/i;Histone H3/j;Histone H3/k;Histone H3/l;HIST2H3A;HIST2H3C;H3F
Swiss	
别名	HIST1H3A/HIST1H3/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B/H3F3C;Histone H3.1/Histone H3.2/Histone H3.3/Histone H3.3 C;H3K27ME2;HIST1H3A;H3FA;HIST1H3B;H3FL;HIST1H3C;H3FC;HIST1H3D;H3FB;HIST1H3E;H3FD;HIST1H3 F;H3FI;HIST1H3G;H3FH;HIST1H3H;H3FK;HIST1H3I;H3FF;HIST1H3J;H3FJ;Histone H3.1;Histone H3/a;Histone H3/b;Histone H3/c;Histone H3/d;Histone H3/f;Histone H3/h;Histone H3/i;Histone H3/j;Histone H3/k;Histone H3/l;HIST2H3A;HIST2H3C;H3F2;H3FM;HIST2H3D;Histone H3.2;Histone H3/m;Histone H3/o;H3F3A;H3.3A;H 3F3;PP781;H3F3B;H3.3B;Histone H3.3;H3F3C;Histone H3.3C;Histone H3.5

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)