

— ABLYBIO, Help Your Research



ASC Rabbit mAb

货号: **AA Y1181**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF IP ELISA
推荐浓度	WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	22kD
实测分子量	22kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.183.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
------	--

基因ID	
基因名	PYCARD;Apoptosis-associated speck-like protein containing a CARD;PYCARD;ASC;CARD5;TMS1;Apoptosis-associated speck-like protein containing a CARD;hASC;Caspase recruitment domain-containing protein 5; PYD and CARD domain-containing protein;Target of methylation-induced silencing 1
Swiss	
别名	PYCARD;Apoptosis-associated speck-like protein containing a CARD;PYCARD;ASC;CARD5;TMS1;Apoptosis-associated speck-like protein containing a CARD;hASC;Caspase recruitment domain-containing protein 5; PYD and CARD domain-containing protein;Target of methylation-induced silencing 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)