

— ABLYBIO, Help Your Research



Trk pan Rabbit mAb

货号: **AA Y1195**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF IP ELISA
推荐浓度	WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	87kD
实测分子量	120-140kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.197.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
------	--

基因ID	
基因名	NTRK1;MTC;TRK;TRKA;NTRK2;TRKB;NTRK3;TRKC;High affinity nerve growth factor receptor;Neurotrophic tyrosine kinase receptor type 1;TRK1-transforming tyrosine kinase protein;Tropomyosin-related kinase A;Tyrosine kinase receptor;Tyrosine kinase receptor A;Trk-A;gp140trk;p140-TrkA;BDNF/NT-3 growth factors receptor;GP145-TrkB;Trk-B;Neurotrophic tyrosine kinase receptor type 2;TrkB tyrosine kinase;Tropomyosin-related kinase B;NT-3 growth factor receptor;GP145-TrkC;Trk-C;Neurotrophic tyrosine kinase rec
Swiss	
别名	NTRK1;MTC;TRK;TRKA;NTRK2;TRKB;NTRK3;TRKC;High affinity nerve growth factor receptor;Neurotrophic tyrosine kinase receptor type 1;TRK1-transforming tyrosine kinase protein;Tropomyosin-related kinase A;Tyrosine kinase receptor;Tyrosine kinase receptor A;Trk-A;gp140trk;p140-TrkA;BDNF/NT-3 growth factors receptor;GP145-TrkB;Trk-B;Neurotrophic tyrosine kinase receptor type 2;TrkB tyrosine kinase;Tropomyosin-related kinase B;NT-3 growth factor receptor;GP145-TrkC;Trk-C;Neurotrophic tyrosine kinase receptor type 3;TrkC tyrosine kinase;;NTRK1;MTC;TRK;TRKA;High affinity nerve growth factor receptor;Neurotrophic tyrosine kinase receptor type 1;TRK1-transforming tyrosine kinase protein;Tropomyosin-related kinase A;Tyrosine kinase receptor;Tyrosine kinase receptor A;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)