

— ABLYBIO, Help Your Research



LRMP Rabbit mAb

货号: **AA1231**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF ICC FC
推荐浓度	WB: 1:1000-5000 IHC: 1:100-300 FC: 1:100-300
理论分子量	
实测分子量	61kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.233.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
------	--

基因ID	
基因名	IRAG2;JAW1;LRMP;Inositol 1,4,5-triphosphate receptor associated 2;Lymphoid-restricted membrane protein;Protein Jaw1; [Cleaved into: Processed inositol 1,4,5-triphosphate receptor associated 2];IRAG2;JAW1;LRMP;Inositol 1,4,5-triphosphate receptor associated 2 ;Lymphoid-restricted membrane protein;Protein Jaw 1; [Cleaved into: Processed inositol 1,4,5-triphosphate receptor associated 2];
Swiss	
别名	IRAG2;JAW1;LRMP;Inositol 1,4,5-triphosphate receptor associated 2;Lymphoid-restricted membrane protein;Protein Jaw1; [Cleaved into: Processed inositol 1,4,5-triphosphate receptor associated 2];IRAG2;JAW1;LRMP;Inositol 1,4,5-triphosphate receptor associated 2 ;Lymphoid-restricted membrane protein;Protein Jaw 1; [Cleaved into: Processed inositol 1,4,5-triphosphate receptor associated 2];

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)