

— ABLYBIO, Help Your Research



PDGFR- α Rabbit mAb

货号: **AA1269**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF IP ELISA
推荐浓度	WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	123kD
实测分子量	190kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.271.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
------	--

基因ID	
基因名	PDGFRA PDGFR2 RHEPDGFRA;Platelet-derived growth factor receptor alpha (PDGF-R-alpha) (PDGFR-alpha) (Alpha platelet-derived growth factor receptor) (Alpha-type platelet-derived growth factor receptor) (CD140 antigen-like family member A) (CD140a antigen) (Platelet-derived growth factor alpha receptor) (Platelet-derived growth factor receptor 2) (PDGFR-2) (CD antigen CD140a);Platelet-derived growth factor receptor alpha;PDGF-R-alpha;PDGFR-alpha;Alpha platelet-derived growth factor receptor;Alpha-t
Swiss	
别名	PDGFRA PDGFR2 RHEPDGFRA;Platelet-derived growth factor receptor alpha (PDGF-R-alpha) (PDGFR-alpha) (Alpha platelet-derived growth factor receptor) (Alpha-type platelet-derived growth factor receptor) (CD140 antigen-like family member A) (CD140a antigen) (Platelet-derived growth factor alpha receptor) (Platelet-derived growth factor receptor 2) (PDGFR-2) (CD antigen CD140a);Platelet-derived growth factor receptor alpha;PDGF-R-alpha;PDGFR-alpha;Alpha platelet-derived growth factor receptor;Alpha-type platelet-derived growth factor receptor;CD140 antigen-like family member A;CD140a antigen;Platelet-derived growth factor alpha receptor;Platelet-derived growth factor receptor 2;PDGFR-2;CD antigen CD140a;

产品验证**实验步骤**

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)