

— ABLYBIO, Help Your Research



PDGFR- β Rabbit mAb

货号: **AA1277**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	IHC: 1:200-1:1000 WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	124kD
实测分子量	190kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.279.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	PDGFRB;Platelet-derived growth factor receptor beta;PDGFRB;PDGFR;PDGFR1;Platelet-derived growth factor receptor beta;PDGF-R-beta;PDGFR-beta;Beta platelet-derived growth factor receptor;Beta-type platelet-derived growth factor receptor;CD140 antigen-like family member B;Platelet-deri
Swiss	
别名	PDGFRB;Platelet-derived growth factor receptor beta;PDGFRB;PDGFR;PDGFR1;Platelet-derived growth factor receptor beta;PDGF-R-beta;PDGFR-beta;Beta platelet-derived growth factor receptor;Beta-type platelet-derived growth factor receptor;CD140 antigen-like family member B;Platelet-deri

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)