

— ABLYBIO, Help Your Research



Rpb1 Rabbit mAb

货号: **AA1288**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA CHIP Cut&Tag
推荐浓度	IHC: 1:50-1:200 WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200 CHIP: 1:50-1:100 Cut&Tag: 1:50-1:100
理论分子量	192kD
实测分子量	250kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.290.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	POLR2A;DNA-directed RNA polymerase II subunit RPB1;POLR2A;POLR2;DNA-directed RNA polymerase II subunit RPB1;RNA polymerase II subunit B1;DNA-directed RNA polymerase II subunit A;DNA-directed RNA polymerase III largest subunit;RNA-directed RNA polymerase II subunit RPB1
Swiss	
别名	POLR2A;DNA-directed RNA polymerase II subunit RPB1;POLR2A;POLR2;DNA-directed RNA polymerase II subunit RPB1;RNA polymerase II subunit B1;DNA-directed RNA polymerase II subunit A;DNA-directed RNA polymerase III largest subunit;RNA-directed RNA polymerase II subunit RPB1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)