

— ABLYBIO, Help Your Research



ERK1/2 Rabbit mAb

货号: **AA1341**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	IHC: 1:2000-1:10000 WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	43kD,41kD
实测分子量	43kD,41kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.343.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	MAPK1/MAPK3;Mitogen-activated protein kinase 3;MAPK3;ERK1;PRKM3;Mitogen-activated protein kinase 3;MAP kinase 3;MAPK 3;ERT2;Extracellular signal-regulated kinase 1;ERK-1;Insulin-stimulated MAP2 kinase;MAP kinase isoform p44;p44-MAPK;Microtubule-associated protein 2 kinase;p
Swiss	
别名	MAPK1/MAPK3;Mitogen-activated protein kinase 3;MAPK3;ERK1;PRKM3;Mitogen-activated protein kinase 3;MAP kinase 3;MAPK 3;ERT2;Extracellular signal-regulated kinase 1;ERK-1;Insulin-stimulated MAP2 kinase;MAP kinase isoform p44;p44-MAPK;Microtubule-associated protein 2 kinase;p

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)