

— ABLYBIO, Help Your Research



# mTOR Rabbit mAb

货号: **AA1343**

## 产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	<b>IHC:</b> 1:200-1:1000 <b>WB:</b> 1:2000-1:10000 <b>IF:</b> 1:200-1:1000 <b>ELISA:</b> 1:5000-1:20000 <b>IP:</b> 1:50-1:200
理论分子量	289kD
实测分子量	260kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.345.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	
基因ID	
基因名	MTOR;Serine/threonine-protein kinase mTOR;MTOR;FRAP;FRAP1;FRAP2;RAFT1;RAPT1;Serine/threonine-protein kinase mTOR;FK506-binding protein 12-rapamycin complex-associated protein 1;FKBP12-rapamycin complex-associated protein;Mammalian target of rapamycin;mTOR;Mechanistic tar
Swiss	
别名	MTOR;Serine/threonine-protein kinase mTOR;MTOR;FRAP;FRAP1;FRAP2;RAFT1;RAPT1;Serine/threonine-protein kinase mTOR;FK506-binding protein 12-rapamycin complex-associated protein 1;FKBP12-rapamycin complex-associated protein;Mammalian target of rapamycin;mTOR;Mechanistic tar

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)