

— ABLYBIO, Help Your Research



# Ubiquitin K48 Rabbit mAb

货号: **AA Y1446**

## 产品信息

反应	Human,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>IHC:</b> 1:200-1:1000 <b>WB:</b> 1:500-1:2000 <b>IF:</b> 1:200-1:1000 <b>ELISA:</b> 1:5000-1:20000
理论分子量	26kD
实测分子量	13-250kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.448.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	
------	--

<b>基因ID</b>	
<b>基因名</b>	UBA52/RPS27A/UBB/UBC;Ubiquitin;UBB;Polyubiquitin-B;UBC;Polyubiquitin-C;RPS27A;UBA80;UBCEP1;Ubiquitin-40S ribosomal protein S27a;Ubiquitin carboxyl extension protein 80;UBA52;UBCEP2;Ubiquitin-60S ribosomal protein L40;CEP52;Ubiquitin A-52 residue ribosomal pr;Ubiquitin K48;
<b>Swiss</b>	
<b>别名</b>	UBA52/RPS27A/UBB/UBC;Ubiquitin;UBB;Polyubiquitin-B;UBC;Polyubiquitin-C;RPS27A;UBA80;UBCEP1;Ubiquitin-40S ribosomal protein S27a;Ubiquitin carboxyl extension protein 80;UBA52;UBCEP2;Ubiquitin-60S ribosomal protein L40;CEP52;Ubiquitin A-52 residue ribosomal pr;Ubiquitin K48;

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)