

— ABLYBIO, Help Your Research



# Transferrin Receptor Rabbit mAb

货号: **AA1480**

## 产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	<b>IHC:</b> 1:200-1:1000 <b>WB:</b> 1:1000-1:5000 <b>IF:</b> 1:200-1:1000 <b>ELISA:</b> 1:5000-1:20000 <b>IP:</b> 1:50-1:200
理论分子量	84kD
实测分子量	84kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.482.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	
基因ID	
基因名	TFRC;Transferrin receptor protein 1 (TR) (TfR) (TfR1) (Trfr) (T9) (p90) (CD antigen CD71) [Cleaved into: Transferrin receptor protein 1, serum form (sTfR)];Transferrin receptor protein 1;TR;TfR;TfR1;Trfr;T9;p90;CD antigen CD71;[Cleaved into: Transferrin receptor protein 1, serum form;sTfR;]
Swiss	
别名	TFRC;Transferrin receptor protein 1 (TR) (TfR) (TfR1) (Trfr) (T9) (p90) (CD antigen CD71) [Cleaved into: Transferrin receptor protein 1, serum form (sTfR)];Transferrin receptor protein 1;TR;TfR;TfR1;Trfr;T9;p90;CD antigen CD71;[Cleaved into: Transferrin receptor protein 1, serum form;sTfR;]

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)