

— ABLYBIO, Help Your Research



Iba1 Rabbit mAb

货号: **AA Y1505**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat, |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC IF IP ELISA |
| 推荐浓度 | IHC: 1:2000-1:4000 WB: 1:1000-1:5000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200 |
| 理论分子量 | 17kD |
| 实测分子量 | 17kD |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.507. |
| 偶联物 | Unconjugated |
| 阳性对照 | |
| 细胞定位 | |
| 纯化 | |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | |
| 基因ID | |
| 基因名 | ;AIF1 G1 IBA1;AIF 1;AIF-1;Aif1;AIF1 protein;AIF1_HUMAN;Allograft inflammatory factor 1;Allograft inflammatory factor 1 splice variant G;allograft inflammatory factor-1 splice variant Hara-1;balloon angioplasty responsive transcription;BART 1;G1;G1 putative splice variant of allograft inflammatory factor 1;IBA 1;IBA1;interferon gamma responsive transcript;Interferon responsive transcript 1;interferon responsive transcript factor 1;ionized calcium binding adapter molecule 1;ionized calcium-binding |
| Swiss | |
| 别名 | ;AIF1 G1 IBA1;AIF 1;AIF-1;Aif1;AIF1 protein;AIF1_HUMAN;Allograft inflammatory factor 1;Allograft inflammatory factor 1 splice variant G;allograft inflammatory factor-1 splice variant Hara-1;balloon angioplasty responsive transcription;BART 1;G1;G1 putative splice variant of allograft inflammatory factor 1;IBA 1;IBA1;interferon gamma responsive transcript;Interferon responsive transcript 1;interferon responsive transcript factor 1;ionized calcium binding adapter molecule 1;ionized calcium-binding adapter molecule 1;ionized calcium-binding adapter molecule;IRT 1;IRT1;Microglia response factor;MRF1;Protein g1; |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)