

— ABLYBIO, Help Your Research



ATF-4 Rabbit mAb

货号: **AA Y1548**

产品信息

反应	Human,Mouse,Rat,
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	IHC: 1:100-1:200 WB: 1:1000-1:5000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	38kD
实测分子量	49kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.550.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	ATF4;Cyclic AMP-dependent transcription factor ATF-4;ATF4;CREB2;TXREB;Cyclic AMP-dependent transcription factor ATF-4;cAMP-dependent transcription factor ATF-4;Activating transcription factor 4;Cyclic AMP-responsive element-binding protein 2;CREB-2;cAMP-responsive element-binding prot
Swiss	
别名	ATF4;Cyclic AMP-dependent transcription factor ATF-4;ATF4;CREB2;TXREB;Cyclic AMP-dependent transcription factor ATF-4;cAMP-dependent transcription factor ATF-4;Activating transcription factor 4;Cyclic AMP-responsive element-binding protein 2;CREB-2;cAMP-responsive element-binding prot

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)