

— ABLYBIO, Help Your Research



TATA Box Binding Protein Rabbit mAb

货号: **AA1558**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat, |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC IF IP ELISA |
| 推荐浓度 | IHC: 1:400-1000 WB: 1:1000-5000 IF: 1:200-1000 ELISA: 1:5000-20000 IP: 1:50-200 |
| 理论分子量 | 38kD |
| 实测分子量 | 38kD |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.560. |
| 偶联物 | Unconjugated |
| 阳性对照 | |
| 细胞定位 | |
| 纯化 | |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | |
| 基因ID | |
| 基因名 | TBP;TATA-box-binding protein (TATA sequence-binding protein) (TATA-binding factor) (TATA-box factor) (Transcription initiation factor TFIID TBP subunit);TBP;GTF2D1;TF2D;TFIID;TATA-box-binding protein;TATA sequence-binding protein;TATA-binding factor;TATA-box factor;Transcription initiation factor TFIID TBP subunit |
| Swiss | |
| 别名 | TBP;TATA-box-binding protein (TATA sequence-binding protein) (TATA-binding factor) (TATA-box factor) (Transcription initiation factor TFIID TBP subunit);TBP;GTF2D1;TF2D;TFIID;TATA-box-binding protein;TATA sequence-binding protein;TATA-binding factor;TATA-box factor;Transcription initiation factor TFIID TBP subunit |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)