

— ABLYBIO, Help Your Research



CDK1/2 (Phospho Thr14) Rabbit mAb

货号: **AA1673**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	IHC: 1:200-1:1000 WB: 1:2000-1:10000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	34kD
实测分子量	34kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.675.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	CDK1/CDK2;Cyclin-dependent kinase 1/2;CDK1;CDC2;CDC28A;CDKN1;P34CDC2;Cyclin-dependent kinase 1;CDK1;Cell division control protein 2 homolog;Cell division protein kinase 1;p34 protein kinase;CDK2;CDKN2;Cyclin-dependent kinase 2;Cell division protein kinase 2;p33 pr
Swiss	
别名	CDK1/CDK2;Cyclin-dependent kinase 1/2;CDK1;CDC2;CDC28A;CDKN1;P34CDC2;Cyclin-dependent kinase 1;CDK1;Cell division control protein 2 homolog;Cell division protein kinase 1;p34 protein kinase;CDK2;CDKN2;Cyclin-dependent kinase 2;Cell division protein kinase 2;p33 pr

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)