

— ABLYBIO, Help Your Research



# TGFβ RI Rabbit mAb

货号: **AA1848**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF IP ELISA
推荐浓度	<b>WB:</b> 1:2000-1:10000 <b>IF:</b> 1:200-1:1000 <b>ELISA:</b> 1:5000-1:20000 <b>IP:</b> 1:50-1:200
理论分子量	56kD
实测分子量	55kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.850.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	
------	--

基因ID	
基因名	TGFBR1 ALK5 SKR4;TGF-beta receptor type-1;TGFBR1;ALK5;SKR4;TGF-beta receptor type-1;TGFR-1;Activin A receptor type II-like protein kinase of 53kD;Activin receptor-like kinase 5;ALK-5;ALK5;Serine/threonine-protein kinase receptor R4;SKR4;TGF-beta type I receptor;Transfor
Swiss	
别名	TGFBR1 ALK5 SKR4;TGF-beta receptor type-1;TGFBR1;ALK5;SKR4;TGF-beta receptor type-1;TGFR-1;Activin A receptor type II-like protein kinase of 53kD;Activin receptor-like kinase 5;ALK-5;ALK5;Serine/threonine-protein kinase receptor R4;SKR4;TGF-beta type I receptor;Transfor

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)