

— ABLYBIO, Help Your Research



MKK3 (Phospho Ser189)/MKK6 (Phospho Ser207) Rabbit mAb

货号: AAY2231

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF IP ELISA
推荐浓度	IHC: 1:100-1:500 WB: 1:1000-1:5000 IF: 1:200-1:1000 ELISA: 1:5000-1:20000 IP: 1:50-1:200
理论分子量	39kD,37kD
实测分子量	39kD,37kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.1233.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	
基因ID	
基因名	MAP2K3;MAP2K6;Serine protease HTRA2 mitochondrial;AW212142;dual specificity mitogen activated protein kinase kinase 3;Dual specificity mitogen activated protein kinase kinase 6;Dual specificity mitogen-activated protein kinase kinase 3;MAP kinase kinase 3;MAP kinase kinase 6;map2k3;MAP2K6;MAPK ERK kinase 3;MAPK/ERK kinase 3;MAPK/ERK kinase 6;MAPKK 3;MAPKK 6;MAPKK3;MAPKK6;MEK 3;MEK 6;MEK3;Mitogen activated protein kinase kinase 3;Mitogen activated protein kinase kinase 6;MKK 3;MKK3;MKK6;mMKK3b;MP
Swiss	
别名	MAP2K3;MAP2K6;Serine protease HTRA2 mitochondrial;AW212142;dual specificity mitogen activated protein kinase kinase 3;Dual specificity mitogen activated protein kinase kinase 6;Dual specificity mitogen-activated protein kinase kinase 3;MAP kinase kinase 3;MAP kinase kinase 6;map2k3;MAP2K6;MAPK ERK kinase 3;MAPK/ERK kinase 3;MAPK/ERK kinase 6;MAPKK 3;MAPKK 6;MAPKK3;MAPKK6;MEK 3;MEK 6;MEK3;Mitogen activated protein kinase kinase 3;Mitogen activated protein kinase kinase 6;MKK 3;MKK3;MKK6;mMKK3b;MP2K3_HUMAN;PRKMK 3;PRKMK3;PRKMK6;Protein kinase, mitogen activated, kinase 6;MAP kinase kinase 6;protein kinase, mitogen-activated, kinase 3;SAPK kinase 2;SAPKK 2;SAPKK 3;SAPKK-2;SAPKK2;SAPKK3;Stress activated protein kinase kinase 2;Stress activated protein kinase kinase 3;Stress-activated protein kinase kinase 2.

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)