

— ABLYBIO, Help Your Research



TGFB3 Antibody

货号: **ABY1003**

产品信息

反应	Human, Mouse, Rat, Monkey
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500-1:3000, IHC 1:50-1:200, IF/ICC 1:100-1:500
理论分子量	47kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	BT-474,K-937,Mouse lung,Mouse skeletal muscle,Mouse heart,Rat lung,Rat heart
细胞定位	[Latency-associated peptide]: Secreted, [Transforming growth factor beta-3]: Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and is involved in embryogenesis and cell differentiation. Defects in this gene are a cause of familial arrhythmogenic right ventricular dysplasia 1.
------	--

基因ID	7043
基因名	TGFB3
Swiss	P10600
别名	ARVD; ARVD1; FLJ16571; LDS5; MGC105479; MGC118722; prepro-transforming growth factor beta-3; RN HF; TGF beta 3; TGF beta3; TGF-beta-3; TGFB 3; Tgfb3; TGFB3_HUMAN; transforming growth factor beta 3; Transforming growth factor beta-3;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)