

— ABLYBIO, Help Your Research



FOXJ1 Antibody

货号: **ABY1031**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500-1:2000, IF/ICC 1:100-1:500, IHC 1:50-1:200
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Nucleus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the forkhead family of transcription factors. Similar genes in zebrafish and mouse have been shown to regulate the transcription of genes that control the production of motile cilia. The mouse ortholog also functions in the determination of left-right asymmetry. Polymorphisms in this gene are associated with systemic lupus erythematosus and allergic rhinitis.
------	---

基因ID	2302
基因名	FOXJ1
Swiss	Q92949
别名	FKHL13; Forkhead box J1; Forkhead box protein J1; Forkhead Homologue 4; Forkhead related protein FKHL 13; Forkhead-related protein FKHL13; Fox J1; Foxj1; FOXJ1_HUMAN; Hepatocyte nuclear factor 3 forkhead homolog 4; HFH 4; HFH-4; HFH4; HNF 3/forkhead homolog 4; MGC35202;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)