

— ABLYBIO, Help Your Research



# CPT1B Antibody

货号: **ABY1049**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500-1:1000, IHC 1:50-1:200, IF/ICC 1:100-1:500
理论分子量	88kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Mitochondrion outer membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene, a member of the carnitine/choline acetyltransferase family, is the rate-controlling enzyme of the long-chain fatty acid beta-oxidation pathway in muscle mitochondria. This enzyme is required for the net transport of long-chain fatty acyl-CoAs from the cytoplasm into the mitochondria. Multiple transcript variants encoding different isoforms have been found for this gene, and read-through transcripts are expressed from the upstream locus that include exons from this gene.
基因ID	1375
基因名	CPT1B
Swiss	Q92523
别名	muscle isoform; Carnitine O palmitoyltransferase 1B; Carnitine O palmitoyltransferase I mitochondrial muscle isoform; Carnitine O palmitoyltransferase I muscle isoform; Carnitine O-palmitoyltransferase 1, muscle isoform; Carnitine O-palmitoyltransferase I; Carnitine palmitoyltransferase 1A (muscle); Carnitine palmitoyltransferase 1B (muscle); Carnitine palmitoyltransferase 1B; Carnitine palmitoyltransferase I like protein; Carnitine palmitoyltransferase I muscle; Carnitine palmitoyltransferase I-like protein; CPT 1B; CPT I; CPT1 M; CPT1 muscle; CPT1-M; Cpt1b; CPT1B_HUMAN; CPT1M; CPTI; CPTI M; CPTI muscle; CPTI-M; CPTIM; FLJ55729; FLJ58750; KIAA1670; M CPT1; M-CPT1; MCCPT1; MCPT1;

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))