

— ABLYBIO, Help Your Research



# CCT6A Antibody

货号: **ABY1068**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500-1:2000,, IHC 1:50-1:200, IF/ICC 1:100-1:500
理论分子量	58kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cytoplasm
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. In addition, several pseudogenes of this gene have been located.
基因ID	908
基因名	CCT6A
Swiss	P40227
别名	Acute morphine dependence related protein 2; Amino acid transport defect complementing; CCT zeta 1; CCT zeta; CCT-zeta; CCT-zeta-1; CCT6; CCTZ; Chaperonin containing T complex subunit 6; Chaperonin containing TCP1 subunit 6A; Chaperonin containing TCP1 subunit 6A (zeta 1); Histidine transport regulator 3; HTR3; MGC126214; MGC126215; MoDP 2; MoDP-2; MoDP2; T complex protein 1 subunit zeta; T complex protein 1 zeta subunit; TCP 1 zeta; TCP-1-zeta; TCP1 zeta; TCP20; TCPZ; TTCP20;

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))