

— ABLYBIO, Help Your Research



FCGR1A Antibody

货号: **ABY1087**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500-1:2000, IHC 1:50-1:200, IF/ICC 1:100-1:500
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	RAW264.7,Mouse liver,Rat lung
细胞定位	Cell membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1.
------	---

基因ID	2209
基因名	FCGR1A
Swiss	P12314
别名	CD64; CD64A; Fc fragment of IgG high affinity Ia receptor (CD64); Fc fragment of IgG, high affinity Ia, receptor for; Fc gamma receptor; Fc gamma receptor I A1; Fc gamma RI; Fc gamma RIA; Fc-gamma RI; Fc-gamma RIA; FCG1; FcgammaRIa; FCGR1; FCGR1_HUMAN; FCGR1A; FcRI; FLJ18345; High affinity immunoglobulin gamma Fc receptor I; IGFR1; IgG Fc receptor I; Immunoglobulin G Fc receptor I; OTTHUMP00000013913;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)