

— ABLYBIO, Help Your Research



# IL21 Antibody

货号: **ABY1088**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500-1:2000,, IHC 1:50-1:200, IF/ICC 1:100-1:500
理论分子量	19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	Raji,Mouse thymus,Rat spleen
细胞定位	Secreted
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the common-gamma chain family of cytokines with immunoregulatory activity. The encoded protein plays a role in both the innate and adaptive immune responses by inducing the differentiation, proliferation and activity of multiple target cells including macrophages, natural killer cells, B cells and cytotoxic T cells. Dysregulation of this gene plays a role in multiple immune-mediated diseases including lupus, psoriasis and chronic inflammatory diseases. Alternatively spliced transcripts encoding multiple isoforms have been observed for this gene.
基因ID	59067
基因名	IL21
Swiss	Q9HBE4
别名	CVID11; IL 21; IL-21; IL21; IL21_HUMAN; Interleukin 21; Interleukin-21; interleukin-21 isoform; Interleukin 21; OTTHUMP00000164088; Za11;

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)