

— ABLYBIO, Help Your Research



CYR61/CCN1 Antibody

货号: **ABY1098**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500-1:2000, IHC 1:50-1:200, IF/ICC 1:100-500
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	HeLa,Mouse heart,Rat lung
细胞定位	Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The secreted protein encoded by this gene is growth factor-inducible and promotes the adhesion of endothelial cells. The encoded protein interacts with several integrins and with heparan sulfate proteoglycan. This protein also plays a role in cell proliferation, differentiation, angiogenesis, apoptosis, and extracellular matrix formation.
------	--

基因ID	3491
基因名	CCN1
Swiss	O00622
别名	3CH61; CCN 1; CCN family member 1; CYR 61; Cyr61; CYR61 protein; CYR61_HUMAN; Cysteine rich angiogenic inducer 61; Cysteine rich heparin binding protein 61; Cysteine-rich angiogenic inducer 61; GIG 1; GI G1; IBP-10; IGF-binding protein 10; Igfbp 10; IGFBP-10; Igfbp10; Insulin like growth factor binding protein 10; Insulin-like growth factor-binding protein 10; Protein CYR61; Protein GIG1;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)