

— ABLYBIO, Help Your Research



Phospho-MLC2 (Ser15) pAb

货号: **ABY1135**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500-1:2000, IF/ICC 1:100-1:500, IHC 1:50-1:200
理论分子量	19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.4.
偶联物	Unconjugated
阳性对照	MCF7,Mouse skeletal muscle
细胞定位	Cytoplasm, myofibril, sarcomere, A band
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Thus gene encodes the regulatory light chain associated with cardiac myosin beta (or slow) heavy chain. Ca ⁺ triggers the phosphorylation of regulatory light chain that in turn triggers contraction. Mutations in t his gene are associated with mid-left ventricular chamber type hypertrophic cardiomyopathy.
------	--

基因ID	4633
基因名	MYL2
Swiss	P10916
别名	Cardiac myosin light chain-2;Cardiac ventricular myosin light chain 2;CMH10;MLC 2v;MLC-2;MLC-2v;MLC 2;MLRV_HUMAN;MYL 2;MYL2;Myosin light chain 2 regulatory cardiac slow;Myosin light polypeptide 2 regulatory cardiac slow;Myosin regulatory light chain 2;Myosin regulatory light chain 2 ventricular/cardiac muscle isoform;Regulatory light chain of myosin;RLC of myosin;Slow cardiac myosin regulatory light chain 2;ventricular/cardiac muscle isoform;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)