

— ABLYBIO, Help Your Research



Phospho-FOXO3A (Ser318) Antibody

货号: **ABY1150**

产品信息

| | |
|-------|---|
| 反应 | Human, Mouse, Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IF/ICC |
| 推荐浓度 | |
| 理论分子量 | 71kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5. |
| 偶联物 | Unconjugated |
| 阳性对照 | 293T |
| 细胞定位 | Cytoplasm, cytosol |
| 纯化 | 亲和纯化 |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|------|--|
| 研究背景 | This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed. |
|------|--|

| | |
|-------|---|
| 基因ID | 2309 |
| 基因名 | FOXO3 |
| Swiss | O43524 |
| 别名 | AF6q21; AF6q21 protein; DKFZp781A0677; FKHR2; FKHL1; FKHL1; FKHL1P2; Forkhead (Drosophila) homolog (rhabdomyosarcoma) like 1; Forkhead box O3; Forkhead box O3A; Forkhead box protein O3; Fork head box protein O3A; Forkhead Drosophila homolog of in rhabdomyosarcoma like 1; Forkhead homolog (rhabdomyosarcoma) like 1; Forkhead in rhabdomyosarcoma like 1; Forkhead in rhabdomyosarcoma-like 1; FOX O3A; FOXO2; foxo3; FOXO3_HUMAN; FOXO3A; MGC12739; MGC31925; |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)