

— ABLYBIO, Help Your Research



Phospho-FOXO1A (Ser256) pAb

货号: **ABY1181**

产品信息

| | |
|-------|---|
| 反应 | Human, Mouse, Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500-1:2000, IHC 1:50-1:200, IF/ICC 1:100-1:500 |
| 理论分子量 | 70kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse liver,Mouse brain,Rat kidney |
| 细胞定位 | Cytoplasm |
| 纯化 | 亲和纯化 |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|------|---|
| 研究背景 | This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma. |
|------|---|

| | |
|-------|---|
| 基因ID | 2308 |
| 基因名 | FOXO1 |
| Swiss | Q12778 |
| 别名 | FKH 1; FKH1; FKHR; Forkhead (Drosophila) homolog 1 (rhabdomyosarcoma); Forkhead box O1; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma; Forkhead, Drosophila, homolog of, in rhabdomyosarcoma; FoxO transcription factor; foxo1; FOXO1_HUMAN; FOXO1A; OTTHUMP0000018301; |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)