

— ABLYBIO, Help Your Research



# Phospho-PINK1 (Ser228) pAb

货号: **ABY1263**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500-1:2000, IHC 1:50-1:200
理论分子量	63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	Rat skeletal muscle
细胞定位	Mitochondrion outer membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from stress-induced mitochondrial dysfunction. Mutations in this gene cause one form of autosomal recessive early-onset Parkinson disease.
------	---

基因ID	65018
基因名	PINK1
Swiss	Q9BXM7
别名	BRPK; FLJ27236; mitochondrial; PARK 6; PARK6; Phosphatase and Tensin Homolog; PINK 1; PINK1; PINK1_HUMAN; Protein kinase BRPK; PTEN induced putative kinase 1; PTEN induced putative kinase protein 1; PTEN-induced putative kinase protein 1; Serine/threonine kinase PINK1 mitochondrial; Serine/threonine protein kinase PINK1 mitochondrial; Serine/threonine-protein kinase PINK1;

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))