

— ABLYBIO, Help Your Research



Phospho-SGK3 (Ser486) pAb

货号: **ABY1272**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500-1:2000, IHC 1:50-1:200
理论分子量	57kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.5.
偶联物	Unconjugated
阳性对照	HL-60,THP-1,Mouse kidney
细胞定位	Cytoplasmic vesicle . Early
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the Ser/Thr protein kinase family and encodes a phosphoprotein with a PX (pho x homology) domain. The protein phosphorylates several target proteins and has a role in neutral amino acid transport and activation of potassium and chloride channels. Alternate transcriptional splice variant s, encoding different isoforms, have been characterized.
------	---

基因ID	23678
基因名	SGK3
Swiss	Q96BR1
别名	CISK; Cytokine independent survival kinase; DKFZp781N0293; Serine/threonine protein kinase Sgk3; Serine/threonine-protein kinase Sgk3; Serum/glucocorticoid regulated kinase 3; Serum/glucocorticoid regulated kinase family member 3; Serum/glucocorticoid regulated kinase like; Serum/glucocorticoid-regulated kinase 3; Serum/glucocorticoid-regulated kinase-like; SGK 2; SGK 3; SGK2; Sgk3; SGK3_HUMAN; SGKL;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)