

Tween 80

货号: **AYB26657**

Tween 80

产品信息

生物活性	Tween 80 (Polysorbate 80), a surfactant, has been widely used as a solvent for pharmacological experiments. Tween 80 can also reduce bacterial attachment and inhibit biofilm formation.
CAS	9005-65-6
中文名称	吐温80
分子量	
体外研究	Tween 80 (Polysorbate 80) is a surfactant which can also reduce bacterial attachment and inhibit biofilm formation. Solubility enhancement as high as 100-fold is achieved using 1% Tween 80 <i>in vitro</i> . The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	Liquid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	
纯度	