

# PEG400

货号: **AYB26660**

---

# PEG400

## 产品信息

生物活性	PEG400 is a strongly hydrophilic polyethylene glycol used as an excellent solvent for a large number of substances. PEG400 is widely used in a variety of pharmaceutical formulations.
CAS	
中文名称	
分子量	
体外研究	
体内研究	
形式	Liquid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	In Vitro: <b>H<sub>2</sub>O : 100 mg/mL</b> (Need ultrasonic) <b>DMSO : 100 mg/mL</b> (Need ultrasonic)
纯度	