

Phosal 50 PG

货号: **AYB26668**

Phosal 50 PG

产品信息

生物活性	Phosal 50 PG is a cosolvent (standardised phosphatidylcholine concentrate). Phosal 50 PG as a carrier for lipophilic compounds, can improve the absorption, efficacy and therapeutic index of the active ingredient.
CAS	774594-96-6
中文名称	
分子量	
体外研究	
体内研究	
形式	Liquid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, protect from light
溶解性	
纯度	