

Nattokinase

货号: **AYB26698**

Nattokinase

产品信息

生物活性	Nattokinase is a potent fibrinolytic enzyme. Nattokinase can break down blood clots by directly hydrolyzing fibrin and plasmin substrate. Nattokinase can be used for the research of cardiovascular diseases.
CAS	133876-92-3
中文名称	
分子量	
体外研究	Nattokinase dissolves 94% blood clot at ten minutes and degrades the fibrin clot after two hours of incubation. The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	<p>In Vitro: H₂O : 100 mg/mL (Need ultrasonic)</p> <p>In Vivo: 请根据您的实验动物和给药方式选择适当的溶解方案。以下溶解方案都请先按照 In Vitro 方式配制澄清的储备液，再依次添加助溶剂： ——为保证实验结果的可靠性，澄清的储备液可以根据储存条件，适当保存；体内实验的工作液，建议您现用现配，当天使用；以下溶剂前显示的百分比是指该溶剂在您配制终溶液中的体积占比；如在配制过程中出现沉淀、析出现象，可以通过加热和/或超声的方式助溶</p> <ul style="list-style-type: none"> 1. 请依序添加每种溶剂： PBS <p>Solubility: 50 mg/mL (Infinity mM); Clear solution; Need ultrasonic *以上所有助溶剂都可在 MCE 网站选购。</p>
纯度	

