


# Polyvinylpyrrolidone

货号: **AYB26700**

## 产品信息

生物活性	Polyvinylpyrrolidone is a compound which has been widely tested and used in human and veterinary medicine as an effective wound healing accelerator and disinfectant when combined with iodine and other compounds.
CAS	9003-39-8
中文名称	聚乙烯吡咯烷酮
分子量	111.1418 (monomer)
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	In Vitro: <b>H<sub>2</sub>O : ≥ 50 mg/mL</b> <b>DMSO : 25 mg/mL</b> (Need ultrasonic) *"≥" means soluble, but saturation unknown.
纯度	≥98.0%