

# Ethoxylated hydrogenated castor oil

货号: **AYB26711**

Ethoxylated hydrogenated castor oil

## 产品信息

生物活性	Ethoxylated hydrogenated castor oil (PEG-40 hydrogenated castor oil) is a combination of synthetic polyethylene glycol (PEG) with natural castor oil. Ethoxylated hydrogenated castor oil can be used to emulsify and solubilize oil-in-water (o/w) emulsions. Ethoxylated hydrogenated castor oil can be used as a cosolvent <i>in vivo</i> .
CAS	61788-85-0
中文名称	聚氧乙烯氢化蓖麻油
分子量	
体外研究	
体内研究	
形式	Oil
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, protect from light
溶解性	In Vitro: <b>H<sub>2</sub>O : ≥ 100 mg/mL</b> *"≥" means soluble, but saturation unknown.
纯度	