

# Potassium dihydrogen phosphate, 99.5%

货号: **AYB26790**

## 产品信息

生物活性	Potassium phosphate monobasic (DiA) is a commonly used in biological assay buffers. Potassium phosphate monobasic is moderate to highly concentrated aqueous solutions of potassium phosphate monobasic for the production of phosphate buffers and other laboratory applications.
CAS	7778-77-0
中文名称	磷酸二氢钾
分子量	136.09
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, sealed storage, away from moisture
溶解性	In Vitro: <b>DMSO : (insoluble or slightly soluble)</b>
纯度	≥99.0%