


# Photo-lysine

货号: **AYB26846**



## 产品信息

生物活性	Photo-lysine, a new lysine-based photo-reactive amino acid, captures proteins that bind lysine post-translational modifications.
CAS	1863117-91-2
中文名称	
分子量	172.19
体外研究	<p>Photo-lysine is designed and synthesized by incorporating a photo-cross-linker (diazirine) into the side chain of natural lysine. Photo-lysine, which is readily incorporated into proteins by native mammalian translation machinery, can be used to capture and identify proteins that recognize lysine post-translational modifications (PTMs), including 'readers' and 'erasers' of histone modifications. Photo-lysine can be incorporated into MDH2 and mediate photo-cross-linking to fix protein-protein interactions in cells. UV irradiation of cells in the presence of photo-lysine induced robust cross-linking of HSP90<math>\beta</math> and HSP60. Photo-lysine has a higher efficiency than photo-leucine for photo-cross-linking of the two chaperone proteins. Photo-lysine enables capture of the heterodimer of proteins Ku70 and Ku80 within a protein complex. Photo-lysine enables identification of histone- and chromatin-binding proteins.</p> <p><b>The accuracy of these methods have not been independently confirmed. They are for reference only.</b></p>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	