

Urokinase (peptidolytic)

货号: **AYB26871**

Urokinase (peptidolytic)

产品信息

生物活性	Urokinase peptidolytic (Urokinase-type plasminogen activator) is a serine protease, an inactive form (zymogen) of the serine protease plasminogen. Activation of plasmin triggers a proteolytic cascade reaction, which in turn participates in thrombolysis or extracellular matrix degradation, implicated in vascular disease and cancer-related research.
CAS	9039-53-6
中文名称	
分子量	
体外研究	Urokinase peptidolytic(uPA) can increase the expression of uPAR protein and mRNA levels in non-small cell lung cancer (NSCLC). The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	