


3'Ome-m7GpppAmpG ammonium

货号: **AYB26877**



产品信息

生物活性	3'Ome-m7GpppAmpG ammonium is a trinucleotide Cap analogue containing a locked nucleic acid (LNA) molecule. 3'Ome-m7GpppAmpG ammonium shows a significant translational efficiency. 3'Ome-m7GpppAmpG ammonium can be used as a potential molecular biology tool in the field of mRNA vaccines and mRNA transfection, such as protein production, gene therapy and anti-cancer immunization.
CAS	
中文名称	
分子量	
体外研究	
体内研究	
形式	Liquid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	-20°C, sealed storage, away from moisture
溶解性	
纯度	≥99.0%