

# Nuclease P1

货号: **AYB26878**

---

## Nuclease P1

### 产品信息

|      |  |
|------|--|
| 生物活性 | Nuclease P1 is a single-stranded specific endonuclease, it hydrolyzes nucleic acids into 5'-mononucleotides and cleaves the single-stranded region of a double-stranded nucleic acid. Nuclease P1 is one of the most well-known single stranded specific nucleases in the field of <b>molecular biology</b> , it is widely used in the pharmaceutical and food industries. |
| CAS  | 54576-84-0   |
| 中文名称 |  |
| 分子量  |  |
| 体外研究 |  |
| 体内研究 |  |
| 形式   | Solid  |
| 运输条件 | Room temperature in continental US; may vary elsewhere.  |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.  |
| 溶解性  |  |
| 纯度   |  |