

Thermolysin

货号: **AYB26894**

Thermolysin

产品信息

生物活性	Thermolysin is a thermostable neutral metalloproteinase enzyme secreted by the Gram-positive bacteria <i>Bacillus thermoproteolyticus</i> . Thermolysin catalyzes the hydrolysis of peptide bonds containing hydrophobic residues.
CAS	9073-78-3
中文名称	嗜热菌蛋白酶
分子量	
体外研究	<p>Thermolysin (100 µg in 1 mL, 64°C for 4 h, in isolated starch granules) digestion results in the selective removal of a group of low-molecular-mass proteins ranging between 10 and 27 kD. Thermolysin can be used for the hydrolysis of bovine liver sarcoplasmic proteins.</p> <p>The accuracy of these methods have not been independently confirmed. They are for reference only.</p>
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	