


DPhPC

货号: **AYB26972**

产品信息

生物活性	DPhPC is a phospholipid used to synthesize bilayer vesicles. DPhPC bilayers do not permit ions to leak in the absence of a pore/ion channel, which can be used for studies on channel proteins.
CAS	207131-40-6
中文名称	
分子量	846.25
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	