


DSPE-PEG-Cyanur

货号: **AYB26991**

产品信息

生物活性	DSPE-PEG-Cyanur is a PEG derivative containing cyanur functional group. DSPE-PEG-Cyanur can be used for PEGylation of protein under mild basic conditions. DSPE-PEG-Cyanur can be used for nanostructured lipid carrier.
CAS	1246304-74-4
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	