


Phosphoric acid (sodium hydrate), $\geq 99.0\%$

货号: **AYB26997**

产品信息

| | |
|------|--|
| 生物活性 | Phosphoric acid (sodium hydrate), $\geq 99.0\%$ is used as a buffer capacity reagent in molecular biology , biochemistry, and chromatography. For the preparation of biological buffers. It is also used to purify antibodies , as a laxative, and in combination with other sodium phosphates |
| CAS | 13472-35-0 |
| 中文名称 | 磷酸二氢钠(二水) |
| 分子量 | 156.01 |
| 体外研究 | |
| 体内研究 | |
| 形式 | Solid |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | 4°C, sealed storage, away from moisture |
| 溶解性 | |
| 纯度 | |