

# Ribonuclease T1

货号: **AYB27049**

## Ribonuclease T<sub>1</sub>

### 产品信息

生物活性	Ribonuclease T1 (Rnase T1), is commonly used in biochemical research. Ribonuclease T1 is an endonuclease that can specifically degrade single stranded RNA. Ribonuclease T1 can form nucleoside 2', 3'-cyclic phosphoric acid intermediates to cut the phosphodiester bond between 3'-guanosine residues and adjacent nucleoside 5'-OH groups to produce 3'-GMP terminal <b>oligonucleotides</b> .
CAS	9026-12-4
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	