


Dioleoylphosphatidylglycerol

货号: **AYB27059**

产品信息

生物活性	Dioleoylphosphatidylglycerol (DOPG) is a naturally occurring phospholipid that can be produced by the combined action of aquaporin 3 (AQP3) and phospholipase D2 (PLD2) . Dioleoylphosphatidylglycerol inhibits macrophage inflammatory mediator production in response to heat shock protein B4 (HSPB4) activation of toll-like receptor 2 (TLR2) .
CAS	62700-69-0
中文名称	
分子量	775.04
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	