

alpha-2,8-Sialyltransferase

货号: **AYB27106**

alpha-2,8-Sialyltransferase

产品信息

生物活性	alpha-2,8-Sialyltransferase (ST8Sia VI; GD3 synthase) is a member of alpha2,8-sialyltransferase (ST8Sia) family, is often used in biochemical studies. alpha-2,8-Sialyltransferase catalyzes elongation of the α 2,8-linked oligo/polysialic acid chain on the Sia residue transferred.
CAS	67339-00-8
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	