


# Bis(tetrabutylammonium) Dichromate

货号: **AYB27135**

## 产品信息

生物活性	Bis(tetrabutylammonium) Dichromate is a biochemical reagent that can be used as a biological material or organic compound for life science related research.
CAS	56660-19-6
中文名称	双(四丁基铵)重铬酸盐 [氧化试剂]
分子量	700.92
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	