


DSPE-PEG

货号: **AYB27166**

产品信息

生物活性	DSPE-PEG is a phospholipids-polymer conjugate that can be used in drug delivery applications. DSPE-PEG is a material for the formulation of nanocarriers for achieving prolonged blood circulation time, improved stability and enhanced encapsulation efficiency.
CAS	145035-96-7
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	