

Glucose 6-phosphate isomerase

货号: **AYB27172**

Glucose 6-phosphate isomerase

产品信息

生物活性	Glucose 6-phosphate isomerase (GPI), i.e., phosphoglucose isomerase/phosphoglucoisomerase (PGI) or phosphohexose isomerase (PHI), is often used in biochemical studies. Glucose 6-phosphate isomerase is a highly conserved glycolytic and dimeric enzyme. Glucose 6-phosphate isomerase catalyzes the reversible isomerization of glucose-6-phosphate to fructose-6-phosphate.
CAS	9001-41-6
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	